

## **User Information Sheet**

locus on quality

| Ronco Item No  | : 951, 961, 971, 981, 991  |                            |  |         |         |     |
|--|--|----------------------------|--|---------|---------|-----|
|  |  |                            |  |         |         |     |
| Brand Name   | RONCO TOUCH (PF)   |                            |  |         |         |     |
| Description<br>Material<br>Donning<br>Characteristics<br>Colour<br>Texture   | <ul> <li>Pink Nitrile Examination Glove</li> <li>Nitrile</li> <li>Powder Free</li> <li>Disposable, Non-Sterile, Ambidextrou</li> <li>Pink (Lavender Rose)</li> <li>Finger textured</li> </ul>  | CHEMO<br>TESTED*           | LANS<br>LANS<br>LANS<br>LANS<br>LANS<br>LANS |         |         |     |
| Specification  | : Size :   | XS                         | S  | М       | L       | XL  |
| opeomodion   | Item Number :  | 951                        | 961  | 971     | <br>981 | 991 |
|  | Length : (mm) $\pm 10$ mm  | 240                        | 240  | 240     | 240     | 240 |
|  | Width : (mm) ± 5 mm  | 75                         | 85   | 95      | 106     | 115 |
|  | Thickness :  |                            |  | 3.5 mil |         |     |
|  | Physical Properties :  |                            |  |         |         |     |
|  | Tensile strength : minimum   | 14 MPa                     |  |         |         |     |
|  | Elongation : minimum   | 500%                       |  |         |         |     |
|  | Pinhole AQL :  | AQL 1.5                    |  |         |         |     |
|  |  |                            |  |         |         |     |
| Packaging  | : C1M = 100 gloves x 10 dispensers = 7   | l,000 gloves               | per case                                     |         |         |     |
| Product Features   | <ul> <li>Extremely flexible for bare hand feel</li> <li>Light weight yet incredibly strong</li> <li>Available in X-Small</li> <li>Pink colour to make the glove visually detectable</li> </ul>   |                            |  |         |         |     |
| Annelisations  |  | <b>A</b>                   |  |         |         |     |
| Applications   | <ul> <li>Dental</li> <li>Examination</li> </ul>  | General Use     Laboratory |  |         |         |     |
|  | Examination     Emergency Medical Services     Aesthetics  |                            |  |         |         |     |
|  | <ul> <li>Long Term Care</li> </ul>   |                            |  |         |         |     |
| Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F ( 37°C ) |  |                            |  |         |         |     |
| Special Notes  | Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s). |                            |  |         |         |     |